



World's fastest autoclave. Class B.

## Another dimension. Same performance.



#### The new definition of sterilization speed.



Time is a relative value. It flows differently under different circumstances. But when you need to have your tools ready, you need them, well, FAST. While other autoclaves might test your patience,

Typical autoclave



Sterilization 134 7 min left

134,0°C STOP 2,10 Bar

## **Enbio Magic Filter: forget about water**



#### Enbio<sup>s</sup>. Speedmaster.

Pouches+

Class B

**Medical Certified** 

Medical device class II b EN-13060 93/42/EEC Directive

Enbio S is compact to fit even a cramped space. Unbelievably short process times will save your procedure in emergency, or allow sterilization of the instruments as you go, between patients.

> Unwrapped FAST / 134°C

7 min

Wrapped Class B / 134°C

15 min

Wrapped Class B / 121°C

30 min







230V/50Hz

2250 W

Unit weight

15 kg

Capacity

Max. packages amount (circa) Max. recommended batch weight 0,5 kg Process chamber volume 2,71

Pressure & temperature

Working pressure Max. pressure Max. process tempreature

2,1 Bar 2,3 Bar 137°C

Dimensions

Width 25 cm Height 16 cm Depth 56 cm Electrical

Power supply Installed power Max. power consumption 10 A Misc

Ingress protection level IP20 Data archiving Water conductivity\*

USB stick < 15 µS/cm

# Over 30,000 autoclaves sold in 6 years.

In just 6 years, Enbio's worldwide community has grown to over 30,000 users, who sterilize instruments in our autoclaves, for the complete safety of their customers and patients.













